

MINUTES



SWODIN __ PACS USER GROUP

DATE: JAN 26 2021

TIME: 1PM-2PM

LOCATION: T/CON (PLEASE REFER TO MEETING INVITATION FOR DETAILS)

Allyse Colafranceschi	x	Kathy Sophocleous	x	Minh Nguyen	x	Yolanda Mundt	
Angela Grundy	x	Kathy Williams		Nick Sauro	x	Trudy DenOuden	
Becki VanderWal	x	Kevin Churilla	x	Sean Milne		Lynn Bender	x
Brenda Murray		Lorelei Thompson		Sheri Bell	x	Caroline Rathwell	
Carolyn MacVicar	x	Luiz Aizawa	x	Spring Robinson	x	Scott King	x
Deryk Howard	x	Mahmoud Badreddine	x	Teresa Crevits	x		
Jun Qin	x	Maureen Hengeveld	x	Wafa Akkila			
Junior Agostini	x	Michelle Harrow	x	Wendi Allinson	x		
Chair: Kathy Sophocleous							
Recorder: Becki Vander Wal							

Topic / Discussion	Action/Decision	Responsible	Due Date
<p>Welcome</p> <ul style="list-style-type: none"> • Approval of Minutes of Nov 26th 2020- Angela Grundy • Review of Action Items from previous meeting- reviewed action items: <ul style="list-style-type: none"> • RA600 follow up – Luiz has completed initial review of options. • Options should be available to sites within the next few months. • Nick S. request options have the ability to import Mammo Tomography 			
<p>2.0 PACS Program Updates  PPT slide for Jan 26.ppt</p> <ul style="list-style-type: none"> • Mid-Line Expansion- <i>Slide 4</i> • Longitudinal Record Project- <i>Slide 4</i> • SWODIN Regional PACS Manual- <i>Slide 4</i> • PACS Expansion to SBGHC <i>Slide 5</i> • Regional Imaging Informatics Specialist- posting <i>Slide 5</i> 			
<p>3.0 PACS Ticket Monitor Update  PPT slide for Jan 26.ppt</p> <ul style="list-style-type: none"> • <i>Slide 6</i> 			
<p>4.0 PACS QC</p> <ul style="list-style-type: none"> • Impact of wrong or deactivate patient ID in PACS <i>Slide 8</i> <ul style="list-style-type: none"> ○ There is a ticket open with Agfa for this issue ○ SWODIN team sends incorrect MRNs off to Agfa weekly for correction ○ These occur when the active MRN in the RIS becomes deactivated in PACS ○ Consequences when they do not match: <ul style="list-style-type: none"> ▪ The overlay in Agfa and ZFP will display the global ID (<i>Slide 8</i>) ▪ URL in Powerchart will not work- send email to London HD and Oneview team will manually correct • PR/SR do not move with merge/fix- <i>Slide 8</i> <ul style="list-style-type: none"> ○ This is generally not an issue, as exams requiring fix /merge are new and do not yet have a PR or SR associated. ○ This is an issue with the migration of old studies that have multiple accession numbers with underscore 1, 2, 3 or 4. 			

Topic / Discussion	Action/Decision	Responsible	Due Date
<ul style="list-style-type: none"> ○ When merging these exams, take note of the number of Images and number of Instances. IF there is a discrepancy between the two numbers and the larger number is on the source exam, do not merge as the PR and/or SR will not move to the target exam. ○ Agfa has this as a fix in a Future Enhancement ● How can you tell physician alias not configured correctly in EI- Slide 8 <ul style="list-style-type: none"> ○ If the physician alias in Agfa does not match the alias in Cerner (which can be distinguished by the prefix of RAD), then the Reading physician name or report author will be inaccurate in the tasks within Agfa. In some cases they will also show as signed off by Administrator, System in the task info. ● Hanging Protocol Logic- Slide 9 <ul style="list-style-type: none"> ○ Nick S- Demo “how to modify a hanging protocol” Slide 9 and 10 ○ Overview of issue: <ul style="list-style-type: none"> ▪ Modalities such as MRI will send series description which can be used to create a specific HP by forcing Agfa HP to hang a given series based on its series description. ▪ The series description must be identical for Agfa HP to select it for the viewport ▪ If the hanging protocol has been created with images from MRI_1, then when MRI_2 is displaying images using this protocol, it may have series that are dithered out as they are not the same as the series name used in the creation of the chosen hanging protocol. Note, it may only be a small space that makes the descriptions different. (slide 9) When there is a discrepancy in the series name, Agfa does not recognize the series and will not hang correctly. ▪ Even within the same modality a tech can rename a series before scanning and any discrepancy will be reflected in the HP not hanging correctly. ▪ ○ Important tips: <ul style="list-style-type: none"> ▪ When in the Hanging Protocol tool, do not use the two buttons at the bottom to modify your HP’s. (Update current HP and Save as New HP). These buttons can corrupt your HP’s. Agfa is looking to disable these buttons in future versions. 	<p>Alias- Follow up with instructions and circle back to SA team</p>	<p>Jun</p>	

Topic / Discussion	Action/Decision	Responsible	Due Date
<ul style="list-style-type: none"> ▪ Do not modify HP's that are system wide, make a personal copy first and then edit ▪ Do- Right click on the HP you want to modify and edit or right click from the clinical side bar and change sort order. ▪ From this menu you can select the series from the dropdown menu ▪ Alternately you can also right click on the HP sequence and edit, this will allow you to drag and drop the series from the clinical sidebar. <ul style="list-style-type: none"> ○ Other points to note with Mammo HP's: <ul style="list-style-type: none"> ▪ Sideways requisitions- this will occurs if the user drags the paper document into a viewport that is purposely inverted for Mammo display. Try dragging to a view port on the second monitor to display in a regular viewport, alternately you can add a HP sequence that has no inherent inversion and drag the paper into that sequence after advancing to it. ▪ Screen caps sent from external sources may be sending without orientation values in dicom tag 0020, 0020, so may also display rotated. Nick S is still investigating this issue with Agfa. Any workarounds or fixes will be shared with the group. ○ Assign context of use for an HP to any User <ul style="list-style-type: none"> ▪ Create HP under your account and give one of the sequences a unique name that you can search it by ▪ Test and then assign to another user in Admin portal ▪ Search Hanging Protocol, find the name and edit ▪ rename the HP from here ▪ assign context of use ▪ call for assistance <ul style="list-style-type: none"> • QC Housekeeping- Slide 11 <ul style="list-style-type: none"> ○ REMINDER – DO NOT USE DELETE ALL FUNCTION – REFER TO PPT 			
<p>5.0 PACS Technical Update</p> <ul style="list-style-type: none"> • Deferred 			

Topic / Discussion	Action/Decision	Responsible	Due Date
6.0 Integrated Systems: <ul style="list-style-type: none"> • Deferred 			
7.0 Round Table: none			

Next Meeting: Mar 30 2021 13:00 – 14:00 Web EX

Legend:  Circulated Material/Attachment provided  Confidential